

SC Tecnoscan probe

Durable scanner for quick pregnancy and backfat detection



SC Tecnoscan probe

S-SWINE — C-CONVEX

Durable scanner for quick pregnancy and backfat detection

Specifications:

- Ultrasound channels: 8
- Battery life: 2.5 hours approx.
- Charging time: 2.5 hours, wireless
- Battery type: Lithium
- Probe: 64 element curved array
- Software: Tecnoscan Connect (Android and iOS)
- Size: 235 x 50 x 70 mm
- Weight: 246 g.
- Waterproof: Yes, splash proof
- Viewing device: Mobile, tablet
- Certifications: CE/FCC/IC/EAC/Anatel

Benefits:

- Image quality: 64 elements
- Durability: electronic probe
- Fully sealed, no wires
- Lightweight ergonomic design
- Innovative low power technology
- Real-time image
- Electronic probe more durable than mechanical ones
- Can connect to any iOS or Android device
- BC scan system
- Backfat and pregnancy scan



App Tecnoscan Connect (Android and iOS)

